

Fitness Gear

Therapeutic Inversion Table



OWNER'S MANUAL

Item #5205

TABLE OF CONTENTS

SERVICE	2
IMPORTANT LABELS	3
IMPORTANT SAFETY INSTRUCTIONS	4
PART DRAWING	5
PART LIST	6
INCLUDED HARDWARE & TOOLS	7
ASSEMBLY	8
SAFETY OPERATING INSTRUCTIONS	17
OPERATION AND ADJUSTMENTS	18
STORAGE	23
WARM UP	24
WARRANTY	25
FAX FORM	26

SERVICE

IMPORTANT: FOR NORTH AMERICA ONLY

**To request product service and order replacement parts, please call our customer service department at:
1-866-924-1688**

Monday through Friday, 8:00 AM-5:00 PM Pacific Standard Time,
or email us at: service@paradigmhw.com
Please visit our website at www.paradigmhw.com.

Please have the following information ready when requesting for service:

- Your name
- Phone number
- Model number
- Serial number
- Part number
- Proof of Purchase

***Before returning this product to the store please contact customer service at the contact number.**

**Paradigm Health & Wellness, Inc.
1189 Jellick Ave, City of Industry, CA 91748, USA**

IMPORTANT LABELS



Paradigm
Model No. 5205
Inversion Table
09/11



CLASSIFIED
UL US
E241212

Inversion Table
UL60601-1
CAN / CSA C22.2 No. 601.1
15HZ

Paradigm Health &
Wellness, Inc.
1189 Jettie Ave.
City of Industry, CA 91748



WARNING
FEET MUST BE SECURED
PROPERLY BEFORE USE

!!WARNING!!
Maximum weight
capacity for this
product is 300LBS/136KGS!!



!!WARNING!!
TO AVOID SERIOUS PHYSICAL INJURY, PLEASE
HOLD ON TO THE HANDLE BAR AND GO BACK
SLOWLY WHILE INVERTING!

FOR ASSISTANCE PLEASE CALL:
1-866-924-1688
08:00 AM- 5:00 PM (PST)
MONDAY THROUGH FRIDAY

WARNING!!
THIS EQUIPMENT IS INTENDED FOR USE BY ADULTS AND CHILDREN 12
YEARS AND OLDER.
POSSIBILITY OF SEVERE INJURY OR DEATH IF USED IMPROPERLY. PLEASE
FOLLOW ALL DIRECTIONS BEFORE USE.
KEEP CHILDREN AWAY!
REPLACE LABELS IF THEY ARE DAMAGED, ILLEGIBLE, OR REMOVED.
PLEASE REVIEW ALL WARNING AND INSTRUCTIONS BEFORE ASSEMBLY
AND USAGE. ALL INSTRUCTIONS SHOULD BE FOLLOWED PRIOR TO
ASSEMBLY AND USAGE.

IMPORTANT SAFETY INSTRUCTIONS


This inversion table was designed and built for optimum safety. However, certain precautions apply whenever you operate the exercise equipment. Be sure to read the entire manual before assembling and operating this equipment. Also, please note the following safety instructions:

1. Consult your physician or other health care professionals before using the inversion table.
2. Always wear proper exercise apparel when using the equipment.
3. If any time you feel faint, light-headed or dizziness while operating the equipment, stop exercise immediately. You should also stop exercising if you are experiencing pain or pressure.
4. Keep children and pets away from the equipment while in use.
5. Only one person should use the equipment at a time.
6. Make sure your equipment is correctly assembled before you use it. Be sure all screws, nuts, and bolts are tightened prior to use.
7. Do not operate this or any exercise equipment if it is damaged.
8. Watch your body: come up slowly, dizziness after a session means you came up too fast. Wait a while after eating before using the inversion table. If you get nauseous, come up as soon as you feel queasy.
9. Always use this equipment on a clear and level surface. Do not use outdoors or near water.
10. Keep hands and feet away from any moving parts. Do not insert any object into any openings.
11. Keep loose clothes, jewelry away from moving parts.
12. **WARNING: ALWAYS HOLD ON TO THE SAFETY HANDLES AND GO BACK SLOWLY WHEN INVERTING. FAILURE TO COMPLY COULD RESULT IN SERIOUS BODILY INJURY.**
13. Children under the age of 12 should not use the following fitness equipment.


WARNING: Before using this equipment you should consult with your personal physician to see if inversion equipment is appropriate for you. Do not use this equipment without your physician's approval. Do not use this equipment if you have any of the following conditions or ailments:

- Extreme obesity
- Glaucoma, retinal detachment or conjunctivitis
- Pregnancy
- Spinal injury, Cerebral Sclerosis, or acutely swollen joints
- Middle ear infection
- High blood pressure, Hypertension, Recent stroke or Transient ischemic attack
- Heart or circulatory disorders for which you are being treated
- Hiatus hernia or Ventral hernia
- Bone weaknesses including Osteoporosis, Unhealed fractures, Modularly pins, or Surgically implanted orthopedic supports
- Use of anti-coagulants including Aspirin in high doses

NOTE: Maximum Weight Capacity for this product is 300lbs/136kgs.




PART DRAWING



PART LIST

No.	Description	Qty	No.	Description	Qty
001	Front U-Frame	1	029	Handlebar	2
002	Rear U-Frame	1	030	Knob	1
003	Adjustable Boom	1	031	Rear Heel Holder	2
004	Bed Frame	1	032	Nylon Strap	1
005	Pivot Arm	2	033	Loop Strap	1
006	Adjustable Instep Frame	1	034	Strap Lock	1
007	Steel Heel Holder Bracket	2	035	Nylon Bed	1
008	Folding Arm	2	036	Foam Grip	2
009	Rod	1	037	Protective Cover	2
010	Bolt M8x23	2	038	Hex Head Bolt M8x23	2
011	Hex Head Bolt M6x47	2	039	Foot Bar	1
012	Phillips Screw M6x30	4	040	Hex Head Bolt M8x50	2
013	Washer Ø20xØ8.5x1.5	16	041	Square End Cap	2
014	Round Plate	1	042	Spring Latch	1
015	Lock Nut M8	8	043	Hex Head Bolt M8x38	2
016	Lock Nut M6	6	044	Rod Cap	2
017	Small Spring Knob	1	045	Plastic Bushing	1
018	Large Spring Knob	1	046	Pad	1
019	Safety Hook	2	047	Pivot Arm Ring	2
020	Rubber Pad	1	048	Screw M6x20	4
021	Oval End Cap	2	049	Right Foot Cap	2
022	Foot Bar End Cap	2	050	Left Foot Cap	2
023	Spring	1	051	Screw M6x25	4
024	Adjustable Instep Frame Square End Cap	1	052	Washer Ø12xØ6.5x1.0	8
025	Round End Cap	2	053	Nut Cap Ø27xØ13.5	2
026	Lower Bed Frame Bushing	2	054	Front Heel Holder	2
027	Washer Ø16xØ6.5x1.0	8	055	VCD	1
028	Upper Bed Frame Bushing	1			


INCLUDED HARDWARE & TOOLS




(11) Hex Head Bolt
M6x47
2 PCS



(13) Washer
Ø20xØ8.5x1.5
12 PCS




(15) Lock Nut M8
6 PCS




(16) Lock Nut M6
2 PCS




(27) Washer
Ø16xØ6.5x1.0
4 PCS



(38) Hex Head Bolt
M8x23
2 PCS




(40) Hex Head Bolt
M8x50
2 PCS




(43) Hex Head Bolt
M8x38
2 PCS



(47) Pivot Arm Ring
2 PCS



(48) Screw M6x20
4 PCS




(51) Screw M6x25
4 PCS




(52) Washer
Ø12xØ6.5x1.0
8 PCS



(53) Nut Cap
Ø27xØ13.5
2 PCS




Phillips Screwdriver
1 PC




Double Open Wrench
M12x14
2 PCS

ASSEMBLY



Tool:



Phillips Screwdriver

Step 1:


Stand up the base of the machine by separating the Front and Rear U-Frames (1, 2). Pull the Front and Rear U-Frames (1, 2) as far apart from each others as possible. Then push down on the middle of the two Folding Arms (8) until they are fully locked down.

Attach the Right & Left Foot Caps (49, 50) to the Front & Rear U-Frame (1, 2) each with one M6x20 Screw (48), one M6x25 Screw (51), and two Ø12xØ6.5x1.0 Washers (52). Tighten screws with the Phillips Screwdriver provided.

Hardware:



(48) Screw M6x20
4 PCS



(51) Screw M6x25
4 PCS




(52) Washer
Ø13xØ6.5x1.0
8 PCS



The product weighs more than 44 lbs and should be assembled and moved by two or more people.


ASSEMBLY



Step 2:



Install two Ø27xØ13.5 Nut Caps (53) onto M8 Lock Nuts (15). Slide one Protective Cover (37) on to each side of the base as shown, and pull down on the Protective Covers (37) until the bottom of the covers are slightly lower than the Folding Arms (8). Use the Velcro straps on the bottom of the Protective Covers (37) to secure the covers to the Folding Arms (8). When the covers are assembled correctly, the Folding Arms (8) should be fully covered by the Protective Covers (37).

ASSEMBLY



Step 3:

Slide the bottom of the Pivot Arms (5) into the brackets that located at each side of the Bed Frame (4), align to the desired hole on the arm with the peg on the bracket. Insert the peg into the hole to lock the Pivot Arms (5) in place. It is recommended that you use the bottom hole on the Pivot Arms (5) until you become more familiar with the equipment.




(47) Pivot Arm Ring
2 PCS

Step 4:

Install the Pivot Arm Rings (47) onto the Pivot Arms (5). Mount the Bed Frame (4) to the Rear U-Frame (2) by inserting the ends of the Pivot Arms (5) into the channels on the plates. The slotted portion of the rollers on the end of the Pivot Arms (5) should be inserted into the channels on the plates.

ASSEMBLY




Step 5:


Slide the Rod (9) through the large round hole on the side of Adjustable Boom (3), and secure the Rod (9) on the Adjustable Boom (3) with one M6x47 Hex Head Bolt (11), one M6 Lock Nut (16), and two Ø16xØ6.5x1.0 Washers (27). Tighten bolt and lock nut with two M12x14 Double Open Wrenches provided. Slide one Steel Heel Holder Bracket (7) and one Rear Heel Holder (31) onto one end of the Rod (9) until the lock tooth is wedged into the slot in the Rod (9) as shown in small figure. Use the same procedure to attach the other Steel Heel Holder Bracket (7) and Rear Heel Holder (31) onto the other end of the Rod (9).

Note: Make sure the lock teeth are wedged into the slots in the Rod (9) to lock the Steel Heel Holder Brackets (7) and Rear Heel Holders (31) in place before use.

Hardware:



(11) Hex Head Bolt
M6x47
1 PC




(16) Lock Nut M6
1 PC



(27) Washer
Ø16xØ6.5x1.0
2 PCS

ASSEMBLY




Step 6:

Slide the Foot Bar (39) into the bottom of the Adjustable Boom (3) and align two of the holes on the Foot Bar (39) with two holes on the boom. Secure the Foot Bar (39) in place using two M8x50 Hex Head Bolts (40), two M8 Lock Nuts (15), and four Ø20xØ8.5x1.5 Washers (13). Tighten bolt and lock nut with two M12x14 Double Open Wrenches provided.


Hardware:



(13) Washer
Ø20xØ8.5x1.5
4 PCS




(15) Lock Nut M8
2 PCS




(40) Hex Head Bolt
M8x50
2 PCS

ASSEMBLY



Tool:




2 Double Open Wrenches
M12x14

Step 7:


Remove the Square End Cap (41) on the back of square bracket of Adjustable Boom (3). Attach the Adjustable Instep Frame (6) to the Adjustable Boom (3) by inserting the Adjustable Instep Frame (6) into the square bracket on the boom. Slide the Adjustable Instep Frame (6) completely into the square bracket, insert the M6x47 Hex Head Bolt (11) with a Ø16xØ6.5x1.0 Washer (27) halfway through the square bracket, slide the M6x47 Hex Head Bolt (11) through the ring at the bottom of the Spring (23), slide the M6x47 Hex Head Bolt (11) through the rest of the square bracket, and secure at the other end with a Ø16xØ6.5x1.0 Washer (27) and M6 Lock Nut (16). Tighten bolt and lock nut with two M12x14 Double Open Wrenches provided. Attach the Square End Cap (41) onto the back of square bracket of Adjustable Boom (3) that was removed.

Note: To slide the Adjustable Instep Frame (6) into the square frame, you must first pull out the Small Spring Knob (17).

Hardware:



(11) Hex Head Bolt
M6x47
1 PC




(16) Lock Nut M6
1 PC



(27) Washer
Ø16xØ6.5x1.0
2 PCS


ASSEMBLY



Step 8:

Pull out the Large Spring Knob (18) and slide the Adjustable Boom (3) into the square bracket on the bottom of the Bed Frame (4) as shown. Slide the Adjustable Boom (3) upward until the desired height on the height scale is just below the bracket on the Bed Frame (4). Lock the Adjustable Boom (3) in place by releasing the Large Spring Knob (18) and sliding the Adjustable Boom (3) up or down slightly until the Large Spring Knob (18) "pops" down into the locked position. For added safety, secured the Knob (30) into the back side of the bracket on the Bed Frame (4) as shown.

ASSEMBLY



Tool:




*2 Double Open Wrenches
M12x14*

Step 9:


Attach the top end of Handlebar (29) onto the Rear U-Frame (2) and Pivot Arm Ring (47) with one M8x23 Hex Head Bolt (38), one M8 Lock Nut (15), and two Ø20xØ8.5x1.5 Washers (13). Tighten bolt and lock nut with two M12x14 Double Open Wrenches provided.

Attach the bottom end of the Handlebar (29) onto the Rear U-Frame (2) with one M8x38 Hex Head Bolt (43), one M8 Lock Nut (15), and two Ø20xØ8.5x1.5 Washers (13). Tighten bolt and lock nut with two M12x14 Double Open Wrenches provided. Use the same procedure to attach the other Handlebar (29) onto the Rear U-Frame (2) and Pivot Arm Ring (47).


Hardware:




(13) Washer
Ø20xØ8.5x1.5
8 PCS



(15) Lock Nut M8
4 PCS




(38) Hex Head Bolt
M8x23
2 PCS




(43) Hex Head Bolt
M8x38
2 PCS

ASSEMBLY



Step 10:

Attach the Nylon Strap (32) to the Strap Lock (34) by inserting the end of the Nylon Strap (32) up through the bottom of the Strap Lock (34), loop the Nylon Strap (32) over the Pre-assembled Loop Strap (33) and down through the Strap Lock (34). Now, loop the strap back over itself, and insert back through the Strap Lock (34), and pull tight to secure. See Diagram.




Step 11:

Attach the Nylon and Loop Straps (32, 33) to the inversion table by hooking the end of the Nylon Strap (32) to the pre-assembled loop on the back of the Bed Frame (4) as shown. Now hook the other end of Loop Strap (33) to the other Pre-assembled loop on the Front U-Frame (1) as shown.


SAFETY OPERATING INSTRUCTIONS

Incorrect



Pivot arm is NOT aligned correctly. The pivot arm is not inserted all the way into the curved slot.



Correct



Make sure the pivot arm is inserted all the way into the slot. Pivot arm is aligned correctly when the groove sits directly on the curved slot and the pivot arm is able to rotate freely.

WARNING: Please make sure both pivot arms are in the same hole to prevent serious injury from occurring.

OPERATION AND ADJUSTMENTS




SHORTEN

LENGTHEN

THE STRAP


For added safety, a nylon strap has been included to restrict the degree of inversion. This strap can be adjusted to different lengths to allow for a greater or lesser degree of inversion. To lengthen the Nylon Strap (32) feed the top end of Nylon Strap (32) into the strap lock, and pull on the lower end of the strap. To shorten the length feed the bottom end of Nylon Strap (32) into the strap lock, and pull on the top end. See Diagram.



ADJUSTING THE BOOM

The Adjustable Boom (3) can be moved to a variety of different positions, in order to accommodate the height of the person on the inversion table. To adjust the boom loosen the knob (30), pull out the Large Spring Knob (18), and slide the boom up or down until the desired height on the height scale is positioned just below the Square Bushing (26). When the boom is in the desired position, simply release the Large Spring Knob (18), slide the boom slightly up or down until the Large Spring Knob (18) locks into place, and tighten the Knob (30).

OPERATION AND ADJUSTMENTS



PIVOT ARMS

The Pivot Arms (5) can be adjusted to allow for a greater or lesser degree of inversion. To adjust the Pivot Arms (5) simply pull out on them until the post is out of the hole, slide them up or down to the desired hole, push in until the post goes through the desired hole. The bottom hole provides the least amount of inversion, while the top hole provides the greatest amount. It is recommended that beginners use the bottom hole until they are familiar with the inversion table.

Note: Both Pivot Arms (5) must be adjusted to the same hole. Trying to adjust the Pivot Arms (5) on two different positions could cause damage to the inversion table, or injury to the user.

THE HANDLEBARS


For added convenience and safety, a set of Handlebars (29) has been added to the inversion table. These Handlebars (29) are located at the top of the Rear U-Frame (2). The Handlebars (29) are there to help you return to the upright position from any degree of inversion. If you wish to return to the upright position, and the bed is moving too slowly, or not moving at all, simply grab the Handlebars (29) and pull on them until you return to the upright position.

Note: The inversion table should always return to the upright position when you move your hands below your waist. If it does not, the inversion table is probably not adjusted correctly to your height.

OPERATION AND ADJUSTMENTS

GENERAL PRECAUTIONS

1. Make sure that the Pivot Arms (5) are locks on the lowest hole for the first few attempts.
2. It is recommended that someone be with you while you are using this inversion table for the first few times.
3. Make sure that the Rear Heel Holders (31) and Front Heel Holders (54) are holding your feet securely.
4. Make sure that the Adjustable Boom (3) is properly set to your height.
5. Make sure that the Adjustable Boom (3) is held securely by both the Large Spring Knob (18) and the Knob (30).
6. Make sure that there is enough room for the bed to rotate completely.



ADJUSTING THE ANGLE HOLDER

1. Pull up on the Small Spring Knob (17), slide the Adjustable Instep Frame (6) completely out of the Adjustable Boom (3).
2. Slide your ankles between the Rear Heel Holders (31) and Front Heel Holders (54) and stand on the foot bar located at the bottom of the Adjustable Boom (3).
3. Pull up on the Small Spring Knob (17), allow the Adjustable Instep Frame (6) to slide back into the Adjustable Boom (3). Push in slightly on the Adjustable Instep Frame (6) until the Rear Heel Holders (31) and Front Heel Holders (54) are around your ankles. Release the Small Spring Knob (17) and adjust the Adjustable Instep Frame (6) slightly until the Small Spring Knob (17) locks into place.
4. Stand upright with your back against the bed and your hands lowered at your sides.


OPERATION AND ADJUSTMENTS

BALANCING THE INVERSION TABLE

The inversion table is like a very sensitively balanced fulcrum. It responds to very slight changes in weight distribution. So, it is very important to make sure that the height is adjusted properly. To do this, mount the inversion table, lock your ankles into the heel holders, and lie back with your hands at your sides. Slowly place your hands across your chest. While in this position, your head should still be above your feet. If your feet are above your head, dismount and adjust the height again.

USING THE INVERSION TABLE

1. Start by lying fully back on the bed with your hands at your side, or resting on your thighs.
2. Keeping your hands close to your body begin to raise your arms slowly allowing the table to rotate backward. Stop, or lower your arms to control the downward rotation of the table.
3. Raise your arms until they are over your head. At this point, the inversion table will be as far back as it can go.
4. As you get more comfortable with the use, rock the bed slowly by moving your arms up and down slowly.
5. It is recommended that the inversion table be used for five or ten minutes each morning, and again each evening.
6. Return to the upright position by slowly moving your hands back down to your thighs.




OPERATION AND ADJUSTMENTS

SUGGESTIONS FOR USE

1. Begin slowly: invert only 15~20 degrees to begin with. Stay inverted only as long as you are comfortable. Return upright slowly.
2. Make gradual changes: increase the angle only if it is comfortable. Increase angle only a few degrees at a time. Increase the time of use 1~2 minutes up to ten over a period of weeks. Add stretching and light exercise only after you are comfortable with inversion.
3. Watch your body. Come up slowly, dizziness after a session means you came up too fast. Wait a while after eating before using table. If you get nauseous, do not fight it, come up as soon as you feel queasy.
4. Keep moving: movement while inverted encourages blood, circulation. Movement may be accomplished by either rhythmic traction or light exercise. Do not exercise strenuously while inverted, limit partial inversion without movement to one or two minutes. Limit full inversion with no movement to only a few seconds.
5. Invert regularly: we recommend two or three times a day depending upon your current condition. Try to schedule it for the same time each day.

STORAGE



For your storage convenience, the inversion table can be folded down to place against a wall, under a bed, or in a storage area.

WARM UP

Quadriceps Stretch

With one hand against a wall for balance, reach behind you and pull your right foot up. Bring your heel as close to your buttocks as possible. Hold for 15 counts and repeat with left foot up.



Inner Thigh Stretch


Sit with the soles of your feet together with your knees pointing outward. Pull your feet as close into your groin as possible.

Gently push your knees towards the floor. Hold for 10 counts.



Toe Touches

Slowly bend forward from your waist, letting your back and shoulders relax as you stretch toward your toes. Reach down as far as you can and hold for 15 counts.



Hamstring Stretches

Sit with your right leg extended. Rest the sole of your left foot against your right inner thigh. Stretch toward your toe as far as possible. Hold for 15 counts. Relax and then repeat with left leg extended.



WARRANTY

The manufacturer warrants to the original purchaser that this product is free from defects in material and workmanship when used for the purpose intended, under the conditions that it has been installed and operated in according to manufacturer's Owner's Manual. Manufacturer's obligation under this warranty is limited to replacing free of charge, any parts which may prove to be defective under normal home use. This warranty does not include any damage caused by improper operation, misuse or commercial application. From the date of purchase, the product is warranted to be free from defects for 1 (one) year. All parts and workmanship, including upholstery, foam, ball bearings, pulleys, cables, shocks, all tension mechanisms, wheels, pedals and hardware are to be free from defects for 90 days. This warranty is offered only to the original owner and is not transferable. Proof of purchase is required. This warranty is offered only to the original owner and is not transferable.

Ordering Replacement Parts

Replacement parts can be ordered by calling or emailing our customer service department

service@paradigmhw.com

1-866-924-1688

Monday through Friday, 8:00 AM - 5:00 PM (PST).

When ordering replacement parts please have the following information ready:

- 1. Owner's Manual**
- 2. Model Number**
- 3. Description of Parts**
- 4. Part Number**
- 5. Date of Purchase**

FAX FORM

PARTS REQUEST FAX FORM

Please fax this form to (1-626-810-2166)

OR YOU CAN EMAIL CUSTOMER SERVICE REQUESTS TO
service@paradigmhw.com

NAME: _____

ADDRESS: _____

CITY _____ STATE _____ ZIP _____

TELEPHONE: (Day) _____

(Night) _____

(Email Address) _____

SERIAL#: _____

MODEL#: _____

PURCHASE DATE: _____

PURCHASE FROM: _____

PART #	DESCRIPTION	QTY

“YOUR ORDER WILL BE PROCESSED WITHIN 3 BUSINESS DAYS”

OFFICIAL USE ONLY

SHIP DATE: _____

TRK #: _____

BACK ORDER: _____