


OWNER'S MANUAL



MASSAGE TABLE




Model # 9096

The specifications of this product may vary from this photo and are subject to change without notice.
IRONMAN, IRONMAN TRIATHLON AND M-DOT are registered trademarks of World Triathlon Corporation.
This product is licensed by the IRONMAN TRIATHLON.

TABLE OF CONTENTS

	<i>Page #</i>
Warning label placement	3
Service	4
Safety precautions	5
Overview drawing	6
Part list	7
Assembly	8
Adjustment	11
Storage	13
Transport the table	13
Maintenance instructions	14
Warranty	15
Part request fax foam	16

WARNING LABEL PLACEMENT



SERVICE

To request for product service and order replacement parts, please call our customer service department at

1-866-924-1688

Monday through Friday, 8:00 am-5:00 pm Pacific Standard Time,
or email at: **service@paradigmhw.com**

Please have the following information ready when requesting for service:

Your name	Phone number	Owner's manual
Model number	Serial number	Part number
Date of Purchase		


***If the product has major defects which prevent it from functioning properly, please return it to the store of purchase within the period allowed by the store.**

Paradigm Health & Wellness, Inc.
1189 Jellick Ave.
City of Industry, CA 91748

SAFETY PRECAUTIONS

- (1) Maintain the table periodically. Make sure all screws are secure and the cable has not deteriorated or frayed. Make sure there are no cracks in the wood structure.
- (2) Do not overload the table beyond the recommended weight limit. This could result in serious injury to the user.
- (3) Do not use alcohol based cleaners when cleaning the vinyl. This could cause damage to the vinyl material. Use warm soapy water and wipe dry.
- (4) Make sure when using the table that all four legs are at the same level. Using the legs at different levels can cause the table to be unstable and cause injury.
- (5) Maximum weight capacity- static weight: 2200 lbs/998 kgs
working weight: 550 lbs/250 kgs

OVERVIEW DRAWING





PART LIST

Part#	Description	Quan.	Part#	Description	Quan.
001	Front Bed with Frame	1	020	Long Cable	2
002	Rear Bed with Frame	1	021	M6x70mm Self Tapping Bolt	4
003	Upper Buckle	2	022	Plastic Washer Ø13xØ6x2t	16
004	Lower Buckle	2	023	Plastic Washer Ø15xØ6x4t	12
005	Handle	2	024	M3x20mm Round Head Self Tapping Screw	10
006	Rubber Foot	4	025	M6 Nylon Nut	4
007	Lower Leg End Cap	4	026	M6x75mm Bolt	4
008	Long Brace	4	027	M6x75mm Self Tapping Bolt	4
009	Cross Brace	2	028	Plastic Washer Ø18xØ6x2t	12
010	Short Brace	4	029	M4x18mm Mod. Truss Head Screw	20
011	Upper Leg	4	030	M6x90mm Bolt	4
012	Lower Leg	4	031	M4x20mm Self Tapping Screw	30
013	Hinge	3	032	U Shape Headrest	1
014	Leg Knob	4	033	Armrest	1
015	End Panel	2	034	Armrest Tube	2
016	M4x20mm Round Head Self Tapping Screw	7	035	1/4"x25.5mm Bolt	4
017	Headrest Bushing Plate	1	036	Headrest Cover	1
018	Face Cushion	1	037	Half Round Bolster	1
019	Short Cable	4			

ASSEMBLY


1

Unlace the Buckles and then open the table so that it is standing on the Rubber Foot side. Take out the Face Cushion, U Shape Headrest, Armrest, and Headrest Cover from underneath of the table.




2

Pull out the left and right pair of Legs and straighten each joint.



3

Lift the table up until all four legs are on the floor.



ASSEMBLY


4

Lift one end of the table and pull out to open table fully.



5

Insert both Right/Left Headrest Sticks with Headrest into the holes of the Headrest Bushing Plate.



6

Place the Face Cushion onto the U Shape Headrest with the designated velcro stickers.

WARNING: HEADREST IS FOR LIGHT RESTING ONLY! PLEASE DO NOT ADD HEAVY PRESSURE WHEN RESTING.



ASSEMBLY


7

Insert the Armrest Tubes into the holes on the Upper Legs.




8

Place the Half Round Bolster on the bed.



P

Place the Headrest Cover on the Face Cushion.



ADJUSTMENT


Adjusting the angle of the Headrest

Pressing up the Locking Handles and then adjust the Headrest to the desired angle. Pressing down the Locking Handles after adjustment.



Adjusting the Table Height

Lay the table on the Rubber Foot side and then remove the Leg Knob.



ADJUSTMENT


Adjusting the Table Height

Remove the Lower Leg and then place the Lower Leg back onto the Upper Leg to the desire hole for adjusting the height of the table.




Tighten the Leg Knob after adjustment.

NOTE: Adjust all of the four Legs to the same height.



STORAGE



1. Lay the table on its side and fold at an angle.
2. Completely fold the Legs flat under the table.
3. Make sure the Cables, U Shape Headrest, Face Cushion, Armrest, and Headrest Cover are inside the table before closing it.

TRANSPORT THE TABLE



The table is easy to carry with the shoulder strap of carry bag.

MAINTENANCE INSTRUCTIONS

You should check your massage table for any kind of wear and tear before each use.

1. Check the legs, beds, cables for wear and tear.
2. Replace damaged and worn components immediately.
3. Keep all damaged equipment out of use until it is repaired.

WARRANTY

Paradigm Inc. warrants to the original purchaser that this product is free from defects in material and workmanship when used for the purpose intended, under the conditions that it has been installed and operated in according to Paradigm's Owner's Manual. Paradigm's obligation under this warranty is limited to replacing free of charge, any parts which may prove to be defective under normal home use. This warranty does not include any damage caused by improper operation, misuse or commercial application. From the date of purchase, the frame is warranted to be free from defects for 1 (one) year. All parts and workmanship, including electronics and its console cases, upholstery, foam, ball bearings, pulleys, cables, shocks, all tension mechanisms, wheels, pedals and hardware are to be free from defects for 90 days. This warranty is offered only to the original owner and is not transferable. Proof of purchase is required.

PART REQUEST FAX FORM

IRONMAN PARTS REQUEST FAX FORM

Please fax this form to (1-626-810-2166)

OR YOU CAN EMAIL CUSTOMER SERVICE REQUESTS TO
service@paradigmhw.com

NAME: _____

ADDRESS: _____

CITY _____ STATE _____ ZIP _____

TELEPHONE: (Day) _____

(Night) _____

(Email Address) _____

SERIAL#: _____

MODEL#: _____

PURCHASE DATE: _____

PURCHASE FROM: _____

PART #	DESCRIPTION	QTY

"YOUR ORDER WILL BE PROCESSED WITHIN 3 BUSINESS DAYS"

OFFICIAL USE ONLY

SHIP DATE: _____

TRK #: _____

BACK ORDER: _____